

## **DreamStation CPAP Pro features**

- · Fixed CPAP therapy modes
- · C-Flex
- Advanced respiratory event detection\*
- · **NEW -** EZ-Start
- · NEW Smart Ramp: 0 to 45 min
- · **NEW** Performance Check
- · Mask fit check
- · Compliance and waveform reports
- · On-board memory
- · NEW remote monitoring solution via Bluetooth®
- · Oximetry module (optional) to monitor therapy efficacy

## Humidification

- · Active humidification controlled by the room environment
- · Pre-heat humidifier capable
- · Heated tube (optional)

## **DreamStation CPAP Pro benefits**

- Daily Progress Feedback displays a simple trend of patients' nightly use to help track progress on the path to adherence.
- **EZ-Start**: automatic, personalized adjustments to CPAP pressure can help patients gradually acclimate to their prescribed level of therapy
- **Connectivity** and remote monitoring options including a standard integrated Bluetooth module
- Performance Check simplifies in-home device evaluation for easy remote troubleshooting
- Patient' data management through detailed compliance and waveform reports

## DreamStation CPAP Pro specifications

Part number	Description
INX400S15	DreamStation CPAP Pro with SD card
INX400H15	DreamStation CPAP Pro with humidifier and SD card

Specifications		
Therapy modes	CPAP, CPAP-Check	
Therapy features	EZ-Start	
Flex pressure relief	0 to 3	
Pressure range	4-20 cm H <sub>2</sub> O	
Smart Ramp time	0 to 45 min (5-min increments)	
Starting ramp pressure	Patient adjustable	
Humidification	Heated humidification: fixed, adaptive, heated tube modes Heated tube (optional)	
Device controls	Color LCD/control wheel/push button	
LCD backlight control	Ambient light sensor	
Filters	Reusable pollen and disposable ultra-fine	
Data storage capacity (mininum)	Display: compliance 1 day, 7-day, 30-day averages SD card > 1 year On-board: up to 6 months	
Altitude compensation	Automatic	
Warranty	2 years	
Dimensions	15,7 (L) $\times$ 19,3 (W) $\times$ 8,4 cm (H) (without humidifier) 29,7 (L) $\times$ 19,3 (W) $\times$ 8,4 cm (H) (with humidifier)	
Weight	1,33 kg (without humidifier), 1,98kg (with humidifier) but including power supply	
Power supply	80W/12V/6,67A	
Electrical requirements	100-240 VAC, 50/60 Hz, 2.0-1.0 A	

Sleep therapy accessories	
Humidification	Part number
DreamStation heated humidifier	EUXH
Water chamber	1122520
Tubing	Part number
Heated tube (15mm)	HT15
Standard tube (15mm)	PR15
Filters	Part number
Pollen (foam) filter, reusable (1/pack)	1122446
Ultra-fine filter, disposable (1/pack)	1122447
Ultra-fine filter, disposable (2/pack)	1122518
Ultra-fine filter, disposable (6/pack)	1122519
Accessories	Part number
DreamStation carrying case	1121162
DreamStation link module	1120293
DreamStation SpO <sub>2</sub> module	1121694
Reuseable finger sensor (clip)	936
Oximetry finger sensor, Adult, Flex	953
Flexiwrap, sens tape, Adult (25/pack)	954A
Reuseable finger sensor (soft), medium	1062050
Reuseable finger sensor (soft), small	1062051
DreamStation, 80W power supply	1118499
Connectivity	Part number
DreamStation cellular/SpO <sub>2</sub> modem	100621C
DreamStation Wifi/SpO <sub>2</sub> modem	100721W





Philips Respironics 1010 Murry Ridge Lane Murrysville, PA 15668 USA

www.philips.com/respironics