

## **Dream**Station

Sleep therapy system



© 2016 Koninklijke Philips N.V. All rights reserved.

Specifications are subject to change without notice. The Bluetooth word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and the use of such marks by Philips Respironics is under license. Wi-Fi is a trademark of Wi-Fi Alliance. www.philips.com/respironics Caution: US federal law restricts these devices to sale by or on the order of a physician.

RRDPGH KW 8/18/16 MCI 4107682 PN 1125341 1010 Murry Ridge Lane, Murrysville, PA 15668 800-345-6443 • 724-387-4000

# Rediscover dreams

# Natural. Comfortable. **Sleep.**

DreamStation supports long-term patient use while creating new efficiencies for your business

### Powerful patient-driven design

DreamStation was designed with feedback from hundreds of patients to provide a simpler PAP experience to help both new and experienced patients use therapy long-term. A suite of intuitive tools can be customized by healthcare professionals to choose the solution that's right for each patient.

#### Acclimating patients in the early days of therapy

**EZ-Start** can help patients gradually adapt to their prescribed level of therapy using an innovative pressure acclimation algorithm with automatic, personalized adjustments to PAP pressure over the first 30 days.\*\*

**SmartRamp** allows users to fall asleep to lower pressures than standard ramp mode. The air pressure remains low until the patient experiences an obstructive respiratory event. Then SmartRamp responds with an auto algorithm to resolve the patient's apneas.\*\*

#### Make progress every day

Each morning patients are greeted with Daily Progress Feedback — a summary of how far they've come to motivate them to stick with their therapy. The intuitive, color display shows a simple trend of their nightly hours of use, followed by a summary of their last 30 days of good nights with more than four hours of use.

### Sleep therapy success\* at your patients' fingertips

Bluetooth® is included in every DreamStation device for easy connectivity to the DreamMapper patient self-management tool. DreamMapper helps patients take an active role in their sleep management by providing daily therapy progress, helpful troubleshooting videos, and goal-setting tools. Patients can stay engaged and get the information they need fast, without having to use your staff's time to resolve common issues.

#### An Encore performance

DreamStation connects to our powerful EncoreAnywhere patient compliance management software. EncoreAnywhere is HIPAA compliant and makes it easy and efficient for sleep professionals to manage patient compliance and therapy.

### Focus your time where it belongs

EncoreAnywhere features Adherence Profiler. It's an optional, intuitive tool that helps you identify patients according to their probability of achieving adherence, helping focus your staff's time on the patients that have the greatest likelihood of success.

\*In a retrospective review conducted by Philips Respironics of approximately 15,000 SystemOne patients, patients who used SleepMapper, which has been rebranded to DreamMapper, demonstrated 22% greater adherence to the therapy than patients who did not use SleepMapper. To see a list of compatible DreamMapper devices, go to www.dreammapper.com/compatible.

\*\*Not available in all DreamStation PAP models and modes.



# Creating efficiencies that make a difference

# Simplify device evaluation

Performance Check simplifies in-home device evaluation for easy remote troubleshooting. Designed to reduce the time and frustration involved in device troubleshooting, this robust tool allows you to guide the patient through an easy remote diagnosis, resulting in a simple "Pass/Fail" result to determine if the device needs to be returned for service.

Performance Check can also generate a report in our EncoreAnywhere patient management software tool, allowing you to troubleshoot patient concerns.

According to our research, 25% of returned devices are found to be operating normally.\* Performance Check is designed to help

reduce the number of normally." Performance Check is designed to help reduce the number of normally functioning devices that are returned for service, saving you time and delivery costs.

# **Reduce** manual adjustments

DreamStation can automatically adjust to changing patient needs, helping reduce costly home visits and letting you focus your attention where it needs to be.

### CPAP-Check

Checks on the user every 30 hours to determine if CPAP pressure is optimal. If not, it automatically adjusts to meet changing patient needs.\*\*

### Auto-Trial

Delivers breath-by-breath auto-CPAP therapy for up to a total 30 days per patient. Transitions to CPAP or CPAP-Check mode at end of Auto-Trial.\*\*

### Opti-Start

Provides optimal starting pressure that can enhance patient comfort and reduce the likelihood of residual events at the beginning of therapy.\*\*

# More cost-effective connectivity options



\*Internal assessment of 2015 competitive CPAP data comparing to ResMed Airsense10/Aircurve10 platform and Fisher & Paykel Icon series platform.

# DreamMapper helps patients take an **active role** in their therapy

With the DreamMapper self-management tool, patients can use the mobile and web-based features to engage with their therapy.



Sync devices Data is easily synced and available through smartphones, tablets, or computer Track progress Patients can monitor the progress of custom goals they set for themselves



Troubleshoot

problems

Troubleshoot common

ney set for themselves



Inspire patients Tools to help encourage adherence Only DreamStation offers Bluetooth<sup>®</sup>, Cellular and Wi-Fi<sup>®</sup> connectivity options. By providing the most comprehensive suite of connectivity choices among leading PAP brands<sup>\*</sup>, you can choose the options that are most cost effective. DreamStation's innovative design has the capability to adapt as technologies continue to advance.

### Making it easy to reach your patients

- $\cdot$  Cellular provides daily data to clinicians and one hour after the patient has turned the device off
- Wi-Fi delivers more data to clinicians more often and provides coverage even in weak cellular areas
- Bluetooth is included in every device for easy connectivity with DreamMapper and Alice NightOne



**View data** Patients can see their therapy data quickly and easily



Set reminders Patients can set reminders for important tasks, including mask, tube and filter cleaning

# The performance and appeal patients want

DreamStation offers a sleek, low-profile design. Its vibrant color screen displays daily progress and features easy-to-use menus.

Dreamstation's top-mounted humidifier hose swivel and front-facing display make it easy for patients to operate while lying down or sitting up in bed. And, it's small and light, making it easy to pack for travel.

30% smaller

The DreamStation PAP is 30% smaller than the Airsense10 PAP.\*\*

**39%** lighter

The DreamStation PAP is 39% lighter than the Airsense10 PAP.\*†



The DreamStation PAP is 63% quieter than ResMed Airsense10.\*\*†

## Hear it from **DreamStation users**

Good-looking piece of equipment! It's modern looking, doesn't look like a medical machine. It's quieter, too. I liked it as soon as I opened it! I liked the way it was laid out and I liked the size of it — it was small, easy to fit on a table along with a lamp and my phone.

\*Philips internal comparison, 2015, data on file. \*\*Based on the average change in sound level during simulated patient breathing, across a range of pressure, expiratory relief, and hu †Does not apply to non-invasive ventilation products



# Insight-driven innovation

sleep diagnostic system



Daily progress feedback Displays previous 3 nights therapy results to track progress

**Detachable humidifier** The humidifier separates to travel even lighter

**Top-mounted hose outlet** Swivel outlet rotates 180° for freedom of movement

Easily share sleep data

with the care team

Simplified design Easy access to water chamber for filling and cleaning

Two-step air Ultra-fine filter

provides additional filtration of very fine particles

# Ordering information

### Sleep therapy

DreamStation Auto CPAP	Part number
DreamStation Auto CPAP	DSX500S11
DreamStation Auto CPAP with humidifier	DSX500H11
DreamStation Auto CPAP with humidifier and heated tube	DSX500T11
DreamStation CPAP Pro	Part number
DreamStation CPAP Pro	DSX400S11
DreamStation CPAP Pro with humidifier	DSX400H11
DreamStation CPAP Pro with humidifier and heated tube	DSX400T11
DreamStation CPAP	Part number
DreamStation CPAP	DSX200S11
DreamStation CPAP with humidifier	DSX200H11
DreamStation CPAP with humidifier and heated tube	DSX200T11
DreamStation Auto BiPAP	Part number
DreamStation Auto BiPAP DreamStation Auto BiPAP	Part number
DreamStation Auto BiPAP	DSX700S11
DreamStation Auto BiPAP DreamStation Auto BiPAP with humidifier DreamStation Auto BiPAP with humidifier	DSX700S11 DSX700H11
DreamStation Auto BiPAP DreamStation Auto BiPAP with humidifier DreamStation Auto BiPAP with humidifier and heated tube	DSX700S11 DSX700H11 DSX700T11
DreamStation Auto BiPAP DreamStation Auto BiPAP with humidifier DreamStation Auto BiPAP with humidifier and heated tube DreamStation BiPAP Pro	DSX700S11 DSX700H11 DSX700T11 Part number
DreamStation Auto BiPAP DreamStation Auto BiPAP with humidifier DreamStation Auto BiPAP with humidifier and heated tube DreamStation BiPAP Pro DreamStation BiPAP Pro	DSX700S11 DSX700H11 DSX700T11 Part number DSX600S11
DreamStation Auto BiPAP DreamStation Auto BiPAP with humidifier DreamStation Auto BiPAP with humidifier and heated tube <b>DreamStation BiPAP Pro</b> DreamStation BiPAP Pro DreamStation BiPAP Pro with humidifier DreamStation BiPAP Pro with humidifier	DSX700S11 DSX700H11 DSX700T11 Part number DSX600S11 DSX600H11
DreamStation Auto BiPAP DreamStation Auto BiPAP with humidifier DreamStation Auto BiPAP with humidifier and heated tube DreamStation BiPAP Pro DreamStation BiPAP Pro DreamStation BiPAP Pro with humidifier DreamStation BiPAP Pro with humidifier and heated tube	DSX700S11 DSX700H11 DSX700T11 <b>Part number</b> DSX600S11 DSX600H11 DSX600T11
DreamStation Auto BiPAP DreamStation Auto BiPAP with humidifier DreamStation Auto BiPAP with humidifier and heated tube DreamStation BiPAP Pro DreamStation BiPAP Pro DreamStation BiPAP Pro with humidifier DreamStation BiPAP Pro with humidifier and heated tube DreamStation BiPAP autoSV	DSX700S11 DSX700H11 DSX700T11 Part number DSX600S11 DSX600H11 DSX600T11 Part number

### **Respiratory** therapy

DreamStation BiPAP AVAPS	Part number
DreamStation BiPAP AVAPS	DSX1130S11
DreamStation BiPAP AVAPS with humidifier	DSX1130H11
DreamStation BiPAP AVAPS with humidifier and heated tube	DSX1130T11
DreamStation BiPAP S/T	Part number
DreamStation BiPAP S/T	Part number DSX1030S11

### Sleep and Respiratory therapy

Humidification	Part number
DreamStation humidifier	DSXH
DreamStation water chamber	1122520
RP-DS Humidifier Dry Box Assembly	1120668
RP-DS Humidifier Flip Lid Seal	1120617
RP-DS Humidifier Dry Box Inlet Seal	1120613
Tubes	Part number
15mm Std Tube – DreamStation - RP	PR15
15mm Heated Tube	HT15
Filters	Part number
Pollen filter, reusable (1 per pack)	1122446
Ultra-fine filter, disposable (1 per pack)	1122447
Ultra-fine filter, disposable (2 per pack)	1122518
Ultra-fine filter, disposable (6 per pack)	1122519

# Sleep and Respiratory therapy

Power Supplies	Part numbe
DreamStation 65W power supply	1120136
DreamStation 80W power supply	1118499
Power cord	1005894
PAP battery kit	Part numbe
PAP Lithium Ion Battery Kit*	1124009
PAP battery kit replacement parts	Part numbe
Li Ion battery, PAP device cable	1121736
Li Ion Battery, Cable Adapter (50 series adapter, #2)	1127799
Li Ion Battery Case	1126947
DreamStation Travel Case	Part numbe
Travel case, 65W power supply w/ USB port	1120135
Travel case, 65W power supply w/ USB port, 5/case	1124301

## Advanced therapies

DreamStation ventilation solutions use clinically proven algorithms that support customized, connected care for patients who require advanced therapy. DreamStation BiPAP autoSV, BiPAP AVAPS, and BiPAP S/T options are designed to offer exceptional comfort, intelligent monitoring, and seamless connectivity between patients and care teams.

Call 1-800-345-6443 or visit www.philips.com/dreamstation for more information.

## Sleep and Respiratory therapy

accessories

Oximetry	Part number
Link module	1120293
Nonin SpO <sub>2</sub> , assembly	1121694
Reuseable finger sensor (clip)	936
Oximetry finger sensor, adult, flex	953
Flexiwrap, sens tape, adult, 25/pk	954A
Reuseable finger sensor (soft), medium	1062050
Reuseable finger sensor (soft), small	1062051
Modems	Part number
Cellular modem, oximetry capable	100600C
Cellular modem, Non-oximetry	100605C
Wi-Fi accessory, oximetry capable	100700W
Wi-Fi accessory, Non-oximetry	100705W
Travel	Part number
Universal PAP/laptop travel briefcase	1114784
DreamStation replacement carrying case	1121162
Shielded DC cord	1120746
Shielded 12V DC power cord system (includes shielded DC cord and battery adapter cable)	1120747
Mask and tubing bedside organizer	Part number
Bedside organizer	1122135
Bedside organizer, 10/case	1124224